

## LS600 Multi-Function Laser Equipment

### Functions:

1. It is effective to remove spots and wrinkles, whiten the skin and shrink pores.
2. Get rid of pigment: Tatoo removal, eyebrows, eyeline, lipline or birth marks.

### Specification:

- |                               |                                     |
|-------------------------------|-------------------------------------|
| - Laser mode:                 | Nd. YAG solid laser                 |
| - Laser wave Length (output): | 1064nm                              |
| - Pulse power:                | ~700mJ                              |
| - Energy density:             | 15mJ/cm <sup>2</sup>                |
| - Pulse Width:                | ~10ns                               |
| - Repeat Frequency:           | 1,2,3,4,5Hz                         |
| - Cooling System:             | Water Cooling in internal recycling |
| - Power:                      | single phase 220V/5A or 110V/10A    |
| - Size:                       | 48*37*36cm                          |
| - N.W:                        | 6Kg                                 |

