

YS-801A Dermabrasion Machine With 9 Heads

The Dermabrasion Machine with 9 heads adopts different thicknesses of sculpture heads, matches with different types of skin and is more convenient used in beauty salons.

The newly developed Miniature Sculpture is a quick, safe and efficient skin renewal technique through shallowly physical polishing with the diamonds of different sizes inbuilt in the sculpture head. The peeling, of which the thickness can be controlled by the beautician along with the texture, has good efficiency for sunshine-aging, pigment deposition and wide pores combined with vacuum suction.

Theory

With the diamond grains in the heads, the Dermabrasion Machine with 9 heads uses vacuum suction to grind the dead skin and dirt out from epidermis, improving skin status, effectively decreasing tiny wrinkles, smoothen some scar like acne and uneven surface of skins. It makes skin more velvet and tender and enhance skins elastic. It is the safest, most effective way of physical skin renewal.

