

PLD-C280 Electromotion Pipe Bed

Specifications

Effective dimension of bed surface:

- Length: 2400mm
- Width: head end: 200mm
- Middle : 400mm foot end:470mm
- Thickness: head end:5mm foot end: 45mm

Movement of Bed Surface:

- Height to the ground: 790~1060mm
- Portrait head side extend: 2000mm max
- Portrait foot side extend: 1300mm max
- Portrait movement: 1600mm
- Transverse movement: 260mm

X-ray absorb quantity - <1.3mmAl

Power Supply:

- Single-phase 220V±22V 50Hz±1Hz
- Capacity: no less than 2kVA

